

Sample Name :KRAS WT Peptide (VVVGAGGVGK)

Cat No. : RP30235

Lot No. : E3800424A

Time Processed :2:15:40

Month-Day-Year Processed :01/08/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: SS-CM-0309

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	8.783	17725	2402	0.296
2	8.967	14700	1727	0.246
3	9.150	36297	5044	0.607
4	9.475	5824994	476805	97.372
5	9.842	47524	6895	0.794
6	20.100	40987	4580	0.685
Total		5982228	497452	100.000