


Mass Spectrum



Sample Information

Month-Day Processed : 06/21/24
Time Processed : 8:46:46
Injection Volume : 0.2
Sample Name : Exendin-4
Cat No. : RP10874
Lot No. : E3845524F
Theoretical MW: 4186.66
Observed MW : 4186.0

Interface

Nebulizing Gas Flow :1.5L/min
CDL Temp :250
Block Temp :200

:ESI

Interface Bias :+4.5 kV
Drying Gas Flow :5 L/min
T.Flow :0.2 ml/min
B.conc :50%H2O/50%MeOH

Equipment

: ZJ21010035