

Sample Name :OVA Peptide (323-339)
Cat No. :RP10610
Lot No. :E3886226G
Time Processed :12:10:21
Month-Day-Year Processed :07/29/2026


Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>
<Detector A>
Column :Inertsil ODS-SP 4.6 x 250 mm
Equipment: ZJ19010324

<Chromatogram>



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	6.842	3815862	550852	95.304
2	7.017	110845	21913	2.768
3	7.292	39396	3445	0.984
4	7.475	21081	3041	0.527
5	7.633	16706	1990	0.417
Total		4003889	581241	100.000