

Sample Name:KRAS G12V Peptide (VVVGAVGVGK)

Cat No. : RP30229

Lot No. : E4068226E

Time Processed :6:48:15

Month-Day-Year Processed :05/16/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ21010026

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	7.683	23930	3492	1.120
2	8.567	2090590	248025	97.854
3	8.900	10402	1621	0.487
4	12.308	11510	1515	0.539
Total		2136431	254653	100.000