

Sample Name :KRAS G12V Peptide (VVVGAVGVGK)

Cat No. : RP30229

Lot No. : E4083124B

Time Processed :2:57:59

Month-Day-Year Processed :03/04/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ19010325

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	9.842	86383	16281	1.219
2	10.116	6740932	1109779	95.140
3	10.594	20716	2853	0.292
4	11.531	46561	3774	0.657
5	12.015	28455	2967	0.402
6	12.365	101920	3737	1.438
7	13.150	9105	1736	0.129
8	15.590	9688	1208	0.137
9	17.704	41502	5818	0.586
Total		7085261	1148155	100.000