


Mass Spectrum



Sample Information

Month-Day Processed : 02/28/24

Time Processed : 9:39:08

Injection Volume : 0.3

Sample Name : TAT peptide

Cat No. : RP20256

Lot No. : E4098224B

Theoretical MW : 1621.91

Observed MW : 1621.6

Interface : ESI

Nebulizing Gas Flow : 1.5L/min

CDL Temp : 250

Block Temp : 200

Equipment

: ZJ21010035

Interface Bias : +4.5 kV

Drying Gas Flow : 5 L/min

T.Flow : 0.2 ml/min

B.conc : 50%H2O/50%MeOH