

Sample Name: Influenza A NP(366-374) Strain A/PR/8/35

Cat No. : RP20267

Lot No. : E4100324B

Time Processed :16:10:49

Month-Day-Year Processed :02/29/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow: 1 ml/min

Wavelength: 220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column : Inertsil ODS-3 4.6 x 250 mm

Equipment: SS-CM-0310

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	8.667	63272	5759	0.694
2	9.033	93681	11610	1.027
3	9.142	13548	3838	0.149
4	10.067	23101	3127	0.253
5	10.217	32837	4219	0.360
6	11.967	8702629	1060050	95.406
7	12.200	99890	16450	1.095
8	12.525	28420	3848	0.312
9	12.733	9353	1790	0.103
10	13.258	8158	1784	0.089
11	13.417	31618	5099	0.347
12	15.067	15163	2435	0.166
Total		9121670	1120008	100.000