

Sample Name :Glucagon-Like Peptide (GLP) II, human

Cat No. :RP10769

Lot No. :E4138926G

Time Processed :23:57:38

Month-Day-Year Processed :07/16/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

| Time  | Module     | Command | Value |
|-------|------------|---------|-------|
| 0.01  | Pumps      | B.Conc  | 5     |
| 25.00 | Pumps      | B.Conc  | 65    |
| 25.01 | Pumps      | B.Conc  | 95    |
| 27.00 | Pumps      | B.Conc  | 95    |
| 27.01 | Pumps      | B.Conc  | 5     |
| 35.00 | Pumps      | B.Conc  | 5     |
| 35.01 | Controller | Stop    |       |

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ21010376

### <Chromatogram>

mV



### <Peak Table>

Detector A Channel 1 220nm

| Peak# | Ret. Time | Area     | Height  | Area%   |
|-------|-----------|----------|---------|---------|
| 1     | 19.608    | 22514    | 4398    | 0.190   |
| 2     | 19.767    | 246409   | 24052   | 2.081   |
| 3     | 20.200    | 107594   | 14062   | 0.909   |
| 4     | 20.517    | 11250974 | 1375489 | 95.006  |
| 5     | 21.000    | 184328   | 18298   | 1.557   |
| 6     | 21.250    | 30541    | 3371    | 0.258   |
| Total |           | 11842359 | 1439669 | 100.000 |