


# Mass Spectrum



|                       |                                       |                       |          |                 |                 |
|-----------------------|---------------------------------------|-----------------------|----------|-----------------|-----------------|
| Sample Information    |                                       | Interface             |          | Equipment       |                 |
| Month-Day Processed : | 07/17/26                              | Nebulizing Gas Flow : | 1.5L/min | Interface Bias  | :ZJ22010150     |
| Time Processed :      | 10:12:21                              | CDL Temp :            | 250      | Drying Gas Flow | : +4.5 kV       |
| Injection Volume :    | 0.6                                   | Block Temp :          | 200      | T.Flow          | :5 L/min        |
| Sample Name :         | Glucagon-Like Peptide (GLP) II, human |                       |          | B.conc          | :0.2 ml/min     |
| Cat No. :             | RP10769                               |                       |          |                 | :50%H2O/50%MeOH |
| Lot No. :             | E4138926G                             |                       |          |                 |                 |
| Theoretical MW :      | 3766.20                               |                       |          |                 |                 |
| Observed MW :         | 3766.4                                |                       |          |                 |                 |