

Sample Name :瑞淑颜™类蛇毒肽
 Cat No. :RP30094CN
 Lot No. :E4260426E
 Time Processed :21:19:34
 Month-Day-Year Processed :05/12/2026
 Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	0
5.00	Pumps	B.Conc	0
20.00	Pumps	B.Conc	30
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	0
35.00	Pumps	B.Conc	0
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm
 Equipment: ZJ21010026

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	6.858	21307	2016	0.243
2	11.308	3423	432	0.039
3	11.617	75698	8763	0.862
4	13.917	8629192	500587	98.249
5	14.983	24652	1175	0.281
6	15.383	28732	1034	0.327
Total		8783004	514005	100.000