


Mass Spectrum



Sample Information

Month-Day Processed :	05/13/26	Equipment	: ZJ21010035
Time Processed :	9:45:19	Interface Bias	: +4.5 kV
Injection Volume :	0.4	Drying Gas Flow	: 5 L/min
Sample Name :	瑞淑颜™ 类蛇毒肽	T.Flow	: 0.2 ml/min
Cat No. :	RP30094CN	B.conc	: 50% H ₂ O / 50% MeOH
Lot No. :	E4260426E		
Theoretical MW :	495.6		
Observed MW :	375.3		