

Sample Name :Hepatitis B Virus Core (128-140)

Cat No. : RP17313

Lot No. : E4319823C

Time Processed :19:08:09

Month-Day-Year Processed :04/01/2023

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	Pump A B.Conc	5
25.00	Pumps	Pump A B.Conc	65
25.01	Pumps	Pump A B.Conc	95
27.00	Pumps	Pump A B.Conc	95
27.01	Pumps	Pump A B.Conc	5
35.00	Pumps	Pump A B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: GR11010440

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	15.396	160209	14271	1.096
2	15.695	14263390	1633436	97.620
3	15.953	91660	11817	0.627
4	16.091	20943	4428	0.143
5	16.798	23395	3305	0.160
6	18.291	51557	6847	0.353
Total		14611156	1674104	100.000