


Mass Spectrum



Sample Information

Month-Day Processed :	04/01/23	Equipment	:ZJ22010150
Time Processed :	23:24:40	Interface Bias	: +4.5 kV
Injection Volume :	0.2	Drying Gas Flow	:5 L/min
Sample Name :	Hepatitis B Virus Core (128-140)	T.Flow	:0.2 ml/min
Sample ID :	U3286WXZG0-1	B.conc	:50%H2O/50%MeOH
Theoretical MW :	1406.63		
Observed MW :	1406.8		