

Sample Name :Glucagon-Like Peptide (GLP) II, human
 Cat No. : RP10769
 Lot No. : E4356224I
 Time Processed :18:36:51
 Month-Day-Year Processed :10/10/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ20010140

<Chromatogram>



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	18.633	186232	9543	1.803
2	19.158	9817743	968179	95.052
3	19.583	236146	25846	2.286
4	19.917	50958	6064	0.493
5	20.433	37714	4507	0.365
Total		10328793	1014139	100.000