

Sample Name :TAT peptide
 Cat No. : RP20256
 Lot No. :E4359025L
 Time Processed :19:44:02
 Month-Day-Year Processed :12/13/2025
 Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	0
5.00	Pumps	B.Conc	0
20.00	Pumps	B.Conc	30
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	0
35.00	Pumps	B.Conc	0
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ24010222

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	10.800	95572	16889	0.433
2	11.158	21417238	2702586	97.040
3	11.475	175988	24658	0.797
4	11.617	112747	11581	0.511
5	11.925	123562	10007	0.560
6	12.383	104737	8173	0.475
7	14.317	40661	4436	0.184
Total		22070504	2778331	100.000