


Mass Spectrum



Sample Information

Month-Day Processed : 04/15/22
 Time Processed : 15:13:18
 Injection Volume : 0.2
 Sample Name : PACAP (1-38), human, ovine, rat
 Cat No. : RP10335
 Lot No. : E4493022C
 Theoretical MW : 4534.26
 Observed MW : 4533.9

Interface

Nebulizing Gas Flow : 1.5L/min
 CDL Temp : 250
 Block Temp : 200

Equipment

Interface Bias : +4.5 kV
 Drying Gas Flow : 5 L/min
 T.Flow : 0.2 ml/min
 B.conc : 50% H₂O/50% MeOH

Equipment

: GK11010007
 : +4.5 kV
 : 5 L/min
 : 0.2 ml/min
 : 50% H₂O/50% MeOH