


Mass Spectrum



Sample Information

Month-Day Processed : 04/01/22
Time Processed : 14:01:03
Injection Volume : 0.1
Sample Name : TAT peptide
Cat No. : RP20256
Lot No. : E4500122C
Theoretical MW : 1621.91
Observed MW : 1622.0

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200
Equipment : ZJ21010035
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH