


Mass Spectrum



Sample Information

Month-Day Processed : 04/04/26

Time Processed : 16:13:23

Injection Volume : 0.3

Sample Name : C-Peptide (3-33), human

Cat No. : RP11116

Lot No. : E4503126D

Theoretical MW: 3020.26

Observed MW : 3019.4

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : ZJ22010150
Interface Bias : -3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH