

Sample Name :TAT peptide
Cat No. :RP20256
Lot No. :E4635825I

Time Processed :16:02:04
Month-Day-Year Processed :09/23/2025

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm
<<LC Time Program>>


Time	Module	Command	Value
0.01	Pumps	B.Conc	0
5.00	Pumps	B.Conc	0
20.00	Pumps	B.Conc	30
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	0
35.00	Pumps	B.Conc	0
35.01	Controller	Stop	

<<Column Performance>>

<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm
Equipment: SS-CM-0310

<Chromatogram>



<Peak Table>

Detector A Channel 1 220nm				
Peak#	Ret. Time	Area	Height	Area%
1	14.333	122998	4550	2.908
2	15.917	2116	562	0.050
3	16.867	10790	1228	0.255
4	18.275	4077713	594382	96.406
5	18.433	16102	9991	0.381
Total		4229719	610714	100.000