


Mass Spectrum



Sample Information

Month-Day Processed : 09/22/25
Time Processed : 3:02:04 PM
Injection Volume : 0.3
Sample Name : TAT peptide
Cat No. : RP20256
Lot No. : E4635825I

Theoretical MW: 1621.91

Observed MW: 1621.6

Interface

Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

:ESI

Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min

Equipment

: ZJ21010035
B.conc : 50%H2O/50%MeOH