


Mass Spectrum



Sample Information

Month-Day Processed : 02/07/23
Time Processed : 22:16:33
Injection Volume : 0.2
Sample Name : TAT peptide
Cat No. : RP20256
Lot No. : E4716923A
Theoretical MW: 1621.91
Observed MW : 1621.6

Interface :ESI
Nebulizing Gas Flow :1.5L/min
CDL Temp :250
Block Temp :200

Equipment
Interface Bias
Drying Gas Flow
T.Flow
B.conc

: ZJ21010035
: +4.5 kV
:5 L/min
:0.2 ml/min
:50%H2O/50%MeOH