

Sample Name :Melittin,honey bee
Cat No. : RP20415
Lot No. : E472223F
Time Processed :2:35:19
Month-Day-Year Processed :06/16/2023
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm
Equipment: SS-CM-0310

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	12.725	12639	1694	0.385
2	12.858	18734	3119	0.571
3	13.050	3173206	440398	96.692
4	13.317	48028	6474	1.463
5	13.400	20667	4300	0.630
6	16.058	8487	1206	0.259
Total		3281761	457191	100.000