


Mass Spectrum



Sample Information

Month-Day Processed : 01/24/24
Time Processed : 15:04:05
Injection Volume : 0.3
Sample Name : DYKDDDDK Peptide
Cat No. : RP10586
Lot No. : E4772024A
Theoretical MW: 1012.97
Observed MW : 1013.0

Interface : ESI
Nebulizing Gas Flow : 1.5 L/min
CDL Temp : 250
Block Temp : 200

Equipment : Z121010035
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH