

Sample Name :Glucagon-Like Peptide (GLP) II, human
 Cat No. : RP10769
 Lot No. : E4865424I

Time Processed :5:20:28

Month-Day-Year Processed :09/23/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ19010141

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	18.629	39565	6000	0.704
2	19.388	132039	31612	2.351
3	19.529	5360833	725558	95.446
4	19.978	49262	2850	0.877
5	20.555	34921	6156	0.622
Total		5616619	772177	100.000