

Sample Name :GRGDSPC
Cat No. :RP20283
Lot No. :E4914726F
Time Processed :18:11:13
Month-Day-Year Processed :06/11/2026
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ17010509

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	4.200	27484	6048	0.343
2	4.917	48086	7095	0.600
3	5.125	23979	2807	0.299
4	5.433	64202	4151	0.801
5	5.783	42236	6531	0.527
6	5.967	7656203	764282	95.466
7	6.550	85901	9280	1.071
8	6.975	71763	13760	0.895
Total		8019854	813953	100.000