

Sample Name :TAT peptide  
Cat No. : RP20256  
Lot No. : E4968622B  
Time Processed :9:06:00 PM  
Month-Day-Year Processed :02/17/2022

Pump A : 0.065% trifluoroacetic in 100% water (v/v)  
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)  
Total Flow:1 ml/min  
Wavelength:220 nm

<<LC Time Program>>

| Time  | Module     | Command       | Value |
|-------|------------|---------------|-------|
| 0.01  | Pumps      | Pump A B.Conc | 5     |
| 25.00 | Pumps      | Pump A B.Conc | 65    |
| 25.01 | Pumps      | Pump A B.Conc | 95    |
| 27.00 | Pumps      | Pump A B.Conc | 95    |
| 27.01 | Pumps      | Pump A B.Conc | 5     |
| 35.00 | Pumps      | Pump A B.Conc | 5     |
| 35.01 | Controller | Stop          |       |


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm  
Equipment: GR11010440

### <Chromatogram>

mV



### <Peak Table>

Detector A Channel 1 220nm

| Peak# | Ret. Time | Area    | Height | Area%   |
|-------|-----------|---------|--------|---------|
| 1     | 8.979     | 1029209 | 83776  | 95.481  |
| 2     | 9.521     | 48714   | 3776   | 4.519   |
| Total |           | 1077922 | 87551  | 100.000 |