


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 03/04/22
Time Processed : 4:34:35 PM
Injection Volume : 0.2
Sample Name : TAT peptide
Cat No. : RP20256
Lot No. : E4968622B
Theoretical MW : 1621.91
Observed MW : 1621.0

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50% H₂O/50% MeOH