


Mass Spectrum



Sample Information

| | | | |
|-------------------------|------------------|-----------------|--------------------|
| Month-Day Processed : | 02/21/23 | Equipment | : ZJ21010035 |
| Time Processed : | 20:12:58 | Interface Bias | : +4.5 kV |
| Injection Volume : | 0.2 | Drying Gas Flow | : 5 L/min |
| Sample Name : | DYKDDDDK Peptide | T.Flow | : 0.2 ml/min |
| Cat No. : | RP10586 | B.conc | : 50% H2O/50% MeOH |
| Lot No. : | E5025723B | | |
| Theoretical MW: 1012.97 | | | |
| Observed MW : 1012.8 | | | |