

Sample Name: Influenza A NP(366-374) Strain A/PR/8/35
Cat No. : RP20267
Lot No. : E5089524E
Time Processed : 12:02:31
Month-Day-Year Processed : 05/16/2024

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow: 1 ml/min
Wavelength: 220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column : Inertsil ODS-SP 4.6 x 250 mm
Equipment: ZJ21010026

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	10.552	15346	2894	0.114
2	11.149	22702	2771	0.169
3	11.525	39482	4062	0.294
4	11.733	22656	2586	0.169
5	11.914	150407	20342	1.121
6	12.177	146458	10791	1.092
7	12.518	35829	4706	0.267
8	12.700	84148	6366	0.627
9	12.955	26832	3126	0.200
10	13.418	12787208	1305349	95.307
11	13.814	85750	6519	0.639
Total		13416817	1369513	100.000