


# Mass Spectrum



## Sample Information

Month-Day Processed : 06/22/26  
Time Processed : 9:58:53  
Injection Volume : 0.3  
Sample Name : TAT peptide  
Cat No. : RP20256  
Lot No. : E5111726F  
Theoretical MW : 1621.91  
Observed MW : 1621.6

Interface : ESI  
Nebulizing Gas Flow : 1.5L/min  
CDL Temp : 250  
Block Temp : 200

Equipment : ZJ21010035  
Interface Bias : +4.5 kV  
Drying Gas Flow : 5 L/min  
T.Flow : 0.2 ml/min  
B.conc : 50%H2O/50%MeOH