


Mass Spectrum



Sample Information		Equipment : ZI21010035	
Month-Day Processed :	04/23/24	Interface Bias :	+4.5 kV
Time Processed :	16:52:43	Drying Gas Flow :	5 L/min
Injection Volume :	0.2	T.Flow :	0.2 ml/min
Sample Name :	DYKDDDDK Peptide	B. conc :	50% H2O/50% MeOH
Cat No. :	RP10586		
Lot No. :	E5121624D		
Theoretical MW:	1012.97		
Observed MW :	1012.8		