

Sample Name :KRAS WT Peptide (VVVGAGGVGK)
Cat No. :RP30235
Lot No. : E5313226E
Time Processed :3:00:21
Month-Day-Year Processed :05/07/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm
Equipment: ZJ20010264

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	8.217	11525	1805	0.447
2	8.408	2489548	341540	96.627
3	8.833	9197	1231	0.357
4	9.042	13588	1565	0.527
5	9.175	12168	1908	0.472
6	9.292	6394	1170	0.248
7	9.483	14197	1606	0.551
8	23.183	19847	1169	0.770
Total		2576465	351995	100.000