


Mass Spectrum



Sample Information

Month-Day Processed :	05/06/26	Interface	:ESI	Equipment	:ZJ22010150
Time Processed :	19:42:09	Nebulizing Gas Flow	:1.5L/min	Interface Bias	: +4.5 kV
Injection Volume :	0.3	CDL Temp	:250	Drying Gas Flow	:5 L/min
Sample Name :	KRAS WT Peptide (VVVGAGGVGK)	Block Temp	:200	T.Flow	:0.2 ml/min
Cat No. :RP30235				B.conc	:50%H2O/50%MeOH
Lot No. : E5313226E					
Theoretical MW : 842.00					
Observed MW : 841.5					