


Mass Spectrum



Sample Information		Interface		Equipment	
Month-Day Processed :	05/06/26	Nebulizing Gas Flow :	1.5L/min	Interface Bias	:ZJ22010150
Time Processed :	18:58:44	CDL Temp :	250	Drying Gas Flow	: +4.5 kV
Injection Volume :	0.3	Block Temp :	200	T.Flow	:5 L/min
Sample Name :	Melittin,honey bee			B.conc	:0.2 ml/min
Cat No. :	RP20415				:50%H2O/50%MeOH
Lot No. :	E5319826E				
Theoretical MW :	2847.50				
Observed MW :	2847.3				