

Sample Name :OVA Peptide (323-339)
Cat No. : RP10610
Lot No. :E5321626E
Time Processed :1:50:06
Month-Day-Year Processed :05/07/2026
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: SS-CM-0310

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	10.658	53081	8996	1.124
2	10.933	56997	8607	1.207
3	11.133	4507875	685783	95.465
4	11.267	85916	14566	1.819
5	11.550	18132	3093	0.384
Total		4722002	721044	100.000