


Mass Spectrum



Sample Information

Acquired by : Gary
Month-Day Processed : 11/07/22
Time Processed : 12:57:12 PM
Injection Volume : 0.6
Sample Name : Jagged-1 (188-204)
Cat No. : RP20331
Lot No. : E5363522J
Theoretical MW : 2107.40
Observed MW : 2107.2

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250|æ
Block Temp : 200|æ

Equipment : Z117010408
Interface Bias : +4.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH