

Sample Name :Glucagon-Like Peptide (GLP) II, human

Cat No. : RP10769

Lot No. : E5372522J

Time Processed :6:11:44 AM

Month-Day-Year Processed :11/17/2022

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ21010375

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	16.771	33161	2251	0.223
2	17.034	15576	1167	0.105
3	18.166	45559	3789	0.306
4	18.708	83490	7288	0.561
5	19.025	232595	21717	1.563
6	19.479	80282	9103	0.540
7	19.571	127433	37733	0.857
8	19.723	14258702	1332349	95.845
Total		14876797	1415397	100.000