


Mass Spectrum



Sample Information

Month-Day Processed :	11/16/22	Equipment	:ZJ22010150
Time Processed :	6:21:09 PM	Interface Bias	: +4.5 kV
Injection Volume :	0.2	Drying Gas Flow	:5 L/min
Sample Name :	Glucagon-Like Peptide (GLP) II, human	T.Flow	:0.2 ml/min
Cat No. :	RP10769	B.conc	:50%H2O/50%MeOH
Lot No. :	E5372522J		
Theoretical MW : 3766.20			
Observed MW : 3766.2			