

Sample Name :OVA Peptide (323-339)

Cat No. :RP10610

Lot No. :E5423926H

Time Processed :7:48:39

Month-Day-Year Processed :08/04/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ19010141

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	8.427	13371566	1404279	96.018
2	8.796	345771	50172	2.483
3	9.095	45143	4388	0.324
4	9.359	61354	6417	0.441
5	9.635	73081	7009	0.525
6	9.896	10683	1602	0.077
7	10.032	18517	1760	0.133
Total		13926115	1475628	100.000