


Mass Spectrum



Sample Information

Month-Day Processed :	07/29/22	Equipment	: ZJ21010035
Time Processed :	11:12:01	Interface Bias	: +4.5 kV
Injection Volume :	0.2	Drying Gas Flow	: 5 L/min
Sample Name :	β-Endorphin, human	T.Flow	: 0.2 ml/min
Cat No. :	RP11344	B.conc	: 50% H ₂ O / 50% MeOH
Lot No. :	E5634122G		
Theoretical MW :	3465.10		
Observed MW :	3464.8		