

Sample Name :Glucagon-Like Peptide (GLP) II, human

Cat No. :RP10769

Lot No. :E5653026C

Time Processed :2:18:10

Month-Day-Year Processed :04/13/2026

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow:1 ml/min

Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	Solvent B Conc.	5
25.00	Pumps	Solvent B Conc.	65
25.01	Pumps	Solvent B Conc.	95
27.00	Pumps	Solvent B Conc.	95
27.01	Pumps	Solvent B Conc.	5
33.00	Pumps	Solvent B Conc.	5
33.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: ZJ19010015

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	17.442	43578	6276	0.249
2	17.708	32606	6125	0.186
3	18.208	44613	5583	0.255
4	18.525	172206	19731	0.984
5	18.842	16708735	1566041	95.476
6	19.267	140429	16093	0.802
7	19.508	358328	31297	2.048
Total		17500494	1651145	100.000