


Mass Spectrum



Sample Information

Month-Day Processed : 11/07/23

Time Processed : 6:55:37 PM

Injection Volume : 0.1

Sample Name : KRAS WT Peptide (VVVGAGGVGK)

Cat No. : RP30235

Lot No. : E5803223K

Theoretical MW : 842.00

Observed MW : 841.4

Interface : ESI
Nebulizing Gas Flow : 1.5L/min
CDL Temp : 250
Block Temp : 200

Equipment : ZJ22010150
Interface Bias : -3.5 kV
Drying Gas Flow : 5 L/min
T.Flow : 0.2 ml/min
B.conc : 50%H2O/50%MeOH