

Sample Name :Melittin,honey bee
Cat No. : RP20415
Lot No. : E5809425B
Time Processed :16:51:31
Month-Day-Year Processed :03/01/2025
Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm

Equipment: ZJ24010222

<Chromatogram>

mV



<Peak Table>

Detector A 220nm

Peak#	Ret. Time	Area	Height	Area%
1	15.233	79387	9597	0.641
2	19.933	22643	2350	0.183
3	20.100	41379	4097	0.334
4	20.483	39117	2956	0.316
5	20.658	65020	6540	0.525
6	20.875	57896	5884	0.468
7	21.250	12045957	1201814	97.337
8	21.717	24152	3733	0.195
Total		12375551	1236972	100.000