

Sample Name :C-Type Natriuretic Peptide (CNP) (1-22), human
Cat No. : RP11110
Lot No. : E5839423F
Time Processed :23:11:20
Month-Day-Year Processed :06/24/2023

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
Total Flow:1 ml/min
Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column :Inertsil ODS-SP 4.6 x 250 mm

Equipment: GK12010012

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	11.818	71040	9788	0.410
2	12.335	132083	18419	0.763
3	12.633	133346	18441	0.770
4	12.829	324141	51330	1.873
5	13.066	16462072	1106195	95.119
6	13.900	131376	11845	0.759
7	14.100	29082	3820	0.168
8	14.616	23625	4162	0.137
Total		17306763	1224001	100.000