


Mass Spectrum



Sample Information

Month-Day Processed :	06/24/23	Interface	:ESI	Equipment	:ZJ22010150
Time Processed :	1:57:28 PM	Nebulizing Gas Flow	:1.5L/min	Interface Bias	: +4.5 kV
Injection Volume :	0.2	CDL Temp	:250	Drying Gas Flow	:5 L/min
Sample Name :	C-Type Natriuretic Peptide (CNP) (1-22), human	Block Temp	:200	T.Flow	:0.2 ml/min
Cat No. :	RP11110	B.conc		B.conc	:50%H2O/50%MeOH
Lot No. :	E5839423F				
Theoretical MW:	2197.60				
Observed MW :	2198.1				