


Mass Spectrum



Sample Information

Month-Day Processed :	02/15/25	Equipment	: ZJ21010035
Time Processed :	13:38:43	Interface Bias	: +4.5 kV
Injection Volume :	0.3	Drying Gas Flow	: 5 L/min
Sample Name :	TAT peptide	T.Flow	: 0.2 ml/min
Cat No. :	RP20256	B.conc	: 50%H2O/50%MeOH
Lot No. :	E5987425B		
Theoretical MW:	1621.91		
Observed MW :	1622.1		