

Sample Name : Influenza A NP(366-374) Strain A/PR/8/35

Cat No. : RP20267

Lot No. : E60439221

Time Processed : 8:32:16

Month-Day-Year Processed : 10/30/2022

Pump A : 0.065% trifluoroacetic in 100% water (v/v)

Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)

Total Flow: 1 ml/min

Wavelength: 220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	Pump A B.Conc	5
25.00	Pumps	Pump A B.Conc	65
25.01	Pumps	Pump A B.Conc	95
27.00	Pumps	Pump A B.Conc	95
27.01	Pumps	Pump A B.Conc	5
35.00	Pumps	Pump A B.Conc	5
35.01	Controller	Stop	

<<Column Performance>>


<Detector A>

Column : Inertsil ODS-SP 4.6 x 250 mm

Equipment: GR11010440

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	9.905	1466612	224633	99.351
2	10.104	9574	2373	0.649
Total		1476186	227006	100.000