


Mass Spectrum



Sample Information

|                       |                                          |                     |           |                 |                              |
|-----------------------|------------------------------------------|---------------------|-----------|-----------------|------------------------------|
| Month-Day Processed : | 10/29/22                                 | Interface           | :ESI      | Equipment       | : ZJ21010035                 |
| Time Processed :      | 17:24:59                                 | Nebulizing Gas Flow | :1.5L/min | Interface Bias  | : +4.5 kV                    |
| Injection Volume :    | 0.2                                      | CDL Temp            | :250      | Drying Gas Flow | :5 L/min                     |
| Sample Name :         | Influenza A NP(366-374) Strain A/PR/8/35 | Block Temp          | :200      | T.Flow          | :0.2 ml/min                  |
| Cat No. :             | RP20267                                  |                     |           | B.conc          | :50%H <sub>2</sub> O/50%MeOH |
| Lot No. :             | E6043922I                                |                     |           |                 |                              |
| Theoretical MW :      | 1026.12                                  |                     |           |                 |                              |
| Observed MW :         | 1026.6                                   |                     |           |                 |                              |