


Mass Spectrum



Sample Information

Month-Day Processed :	05/31/22	Equipment	: GK11010007
Time Processed :	20:30:10	Interface Bias	: +4.5 kV
Injection Volume :	0.2	Drying Gas Flow	:5 L/min
Sample Name :	TAT peptide	T.Flow	:0.2 ml/min
Cat No. :	RP20256	B.conc	:50%H2O/50%MeOH
Lot No. :	E6169322E		
Theoretical MW:1621.91			
Observed MW :1621.6			