

Sample Name :Hepatitis B Virus Core (128-140)
 Cat No. :RP17313
 Lot No. :E6253222J
 Time Processed :13:53:56
 Month-Day-Year Processed :10/31/2022

Pump A : 0.065% trifluoroacetic in 100% water (v/v)
 Pump B : 0.05% trifluoroacetic in 100% acetonitrile (v/v)
 Total Flow:1 ml/min
 Wavelength:220 nm

<<LC Time Program>>

Time	Module	Command	Value
0.01	Pumps	B.Conc	5
25.00	Pumps	B.Conc	65
25.01	Pumps	B.Conc	95
27.00	Pumps	B.Conc	95
27.01	Pumps	B.Conc	5
35.00	Pumps	B.Conc	5
35.01	Controller	Stop	


<<Column Performance>>

<Detector A>

Column :Inertsil ODS-3 4.6 x 250 mm
 Equipment: ZJ19010141

<Chromatogram>

mV



<Peak Table>

Detector A Channel 1 220nm

Peak#	Ret. Time	Area	Height	Area%
1	10.260	23347	1917	0.152
2	10.968	33407	2482	0.218
3	11.191	61933	4414	0.404
4	11.904	257451	15792	1.680
5	12.464	14860020	1082838	96.956
6	13.151	33041	2499	0.216
7	13.711	36215	1879	0.236
8	23.854	21207	737	0.138
Total		15326622	1112558	100.000